

Airplane Wheelchair Entry

Chair Force One
December 2nd, 2025



Stanford | Mechanical Engineering
Senior Capstone Program

Project Members:



Vicky

Haley

Brad

Rodolfo

Team Advisors & Support



Coach:
Bill




CA:
Tori



Liaison:
Sean Auer

Project Motivation:



70% of wheelchair users said they forgo air travel due to accessibility concerns.



Project Objective:

Design a **cabin layout specific** solution
that allows **manual wheelchairs** users to **independently board and exit**
an aircraft, move through its aisle, and secure themselves
on **standard aircraft seat**.



Who are our users?

Primary user:

- A traveling manual wheelchair user
- Has upper body functionality
- Weight under 300 lbs



Secondary users:

- Airline and airport staff
- Primary user's travel companions



High Priority Requirements

UR-1-1	Essential actions can be completed by primary user	ER-1-1	Actions require no more than 12 lbf/hand
UR-1-2	Narrow enough to fit in B737 aisle	ER-1-2	In the aisle , wheelchair width < 18 inches
UR-1-5	Solution must be stowable in the aircraft cabin	ER-1-5	Must fit into a space that is 54 x 22 x 14 inches
UR-3-1	Wheelchair withstands static loading conditions	ER-3-1a	300 lb weight at seat without yielding/fracture
		ER-3-1b	Backrest withstands 105 lbf at top point
		ER-3-1c	Solution must not tip with a max of 15 degrees tilt in all directions
UR-4-1	Power system must conform to FAA standards	ER-4-1	Energy capacity ≤ 300 Wh for Li-Ion Battery

High Priority Requirements

UR-1-1	Essential actions can be completed by primary user	ER-1-1	Actions require no more than 12 lbf/hand
UR-1-2	Narrow enough to fit in B737 aisle	ER-1-2	In the aisle , wheelchair width < 18 inches
UR-1-5	Solution must be stowable in the aircraft cabin	ER-1-5	Must fit into a space that is 54 x 22 x 14 inches
UR-3-1	Wheelchair withstands static loading conditions	ER-3-1a	No yielding/fracture with 300 lb weight on seat
		ER-3-1b	Backrest withstands 105 lbf at top point
		ER-3-1c	Solution must not tip with a max of 15 degrees tilt in all directions
UR-4-1	Power system must conform to FAA standards	ER-4-1	Energy capacity ≤ 300 Wh for Li-Ion Battery

High Priority Requirements

UR-1-1	Essential actions can be completed by primary user	ER-1-1	Actions require no more than 12 lbf/hand
UR-1-2	Narrow enough to fit in B737 aisle	ER-1-2	In the aisle , wheelchair width < 18 inches
UR-1-5	Solution must be stowable in the aircraft cabin	ER-1-5	Must fit into a space that is 54 x 22 x 14 inches
UR-3-1	Wheelchair withstands static loading conditions	ER-3-1a	No yielding/fracture with 300 lb weight on seat
		ER-3-1b	Backrest withstands 105 lbf at top point
		ER-3-1c	Solution must not tip with a max of 15 degrees tilt in all directions
UR-4-1	Power system must conform to FAA standards	ER-4-1	Energy capacity ≤ 300 Wh for Li-Ion Battery


The Foldable, Dual Mode Wheelchair

An adaptation to a COTS foldable manual wheelchair

Overview of 3 Concepts from DR7

Problem: How to get a wheelchair user from the jet bridge to their aircraft seat, and back.

1






2




3





Concept 1: The Electric Air Chair



Concept 2: Manual Wheelchair with Solo Struts



Concept 3: Manual Wheelchair with Telescoping Struts




DR7 Concepts



DR7 Concepts

Applicable ERs	Weight	Electric AirChair	Manual WC w/ Solo Struts	Manual WC w/ Telescoping Strut
ER-1-1a: Propeling takes no more than 12 lbf per hand	3	3	2	2
ER-1-2: Wheelchair width in aisle < 18 inches	3	3	3	3
ER-1-5: Wheelchair in stowed configuration fits in 54 x 22 x 14 inch volume	3	1	3	2
ER-1-6: All actuations take < 15 seconds to deploy	1	3	2	2
ER-3-1a: Wheelchair must support 300 lb weight limit without yield/failure	3	1	2	3
ER-3-1c: Tilt 15 degrees without tip-over with 170 lb weight in seat	3	2	3	3
ER-4-1d: Wheelchair weighs < 50 lbs	2	3	2	1
ER-4-2a: Battery must be < 300 Wh capacity	3	1	3	3
	Total	42	54	52



Solo Struts



DR7 Feedback



A Foldable,
Dual Mode
Wheelchair



The Linear Actuators

- Relevant ERs

- ER-3-1a: Integrated components cannot yield with 300lbs
- ER-3-1c: <15 degree tilt on wheelchair at all times
- ER-4-2a: Lithium ion batteries < 300Wh

- Main Risks:

- **Risk of tipping (loss of synchronization)**



Feedback - Yes

Standard Feedback

Voltage - 12 VDC


12 VDC 24 VDC

Stroke - 3 inch



1 inch 2 inch 3 inch 4 inch 6 inch 8 inch 12 inch
24 inch 40 inch

Force - 169 lbs

16 lbs 56 lbs 169 lbs




Attachment of Linear Actuator to Wheelchair





- Relevant ERs: ER 3-1a & 3-1c
- Main Risk: ***Twisting about connection point***

How are we going to attach the actuator to the wheel?

→ Attachment with a pin hole





→ ER 3-1a: 300 lb weight



Fitting Wheelchair in the Overhead Space

- High Priority ERs
 - ER-1-5: Stowed wheelchair fits in 54 x 22 x 14 inch volume
 - ER-3-1b: Backrest withstands 105 lbf of horizontal aft forces
- Original wheelchair, ER-1-5 unsatisfied → split backrest
- Main Risk: ***The removable backrest connection***



Stress Analysis + FEA of Backrest


Connector Only

$$M = Fd = (52.5 \text{ lbf})(7 \text{ in}) = 367.5 \text{ lb-in}$$

$$I = \frac{\pi}{64} d_o^4 = 0.01655 \text{ in}^4$$


$$\sigma_{max} = \frac{Md}{2I} \approx 8460 \text{ psi}$$

$$\sigma_{yield, 6061} \approx 40,000 \text{ psi} > \sigma_{max}$$



6061 Al \rightarrow $FS_{MIN} > 4$

Backrest Tubes



*Unrealistically high stresses at bolt contacts, simulation refinement required

AISI 1060 Steel \rightarrow $FS_{MIN} \approx 1.7$
(bulk tube stress)

Ethics of Our Design

Consideration

Safety for primary and secondary users (ASME Canon 1)

- {
- Sufficient safety factors (Medium ER)
 - Design across multiple failure modes
 - Wheelchair max weight = 50 lbs
 - Reduce physical strain on secondary users




Ethics of Our Design

Consideration

Balancing client interests of Boeing, airlines, and users (ASME Canon 4)

- Comfort & independence vs. Efficiency (time & cost)




Immediate Next Steps

Subassembly	Action Item
Linear Actuators	<ul style="list-style-type: none">• Decide between rechargeable vs. single-use battery• Purchase hardware: linear actuators, microcontroller, battery
Linear Actuators & WC Frame Attachment	<ul style="list-style-type: none">• Redesign to reduce risk of twisting about frame connection point
Linear Actuator & Wheel Attachment	<ul style="list-style-type: none">• Finalize material selection (CF vs. PLA)
Removable Backrest	<ul style="list-style-type: none">• Refine FEA simulations to improve contacts• Finalize material selection (PC-CF vs. Al)• Develop a full backrest testing procedure

Gantt Chart

1			
2			
3			
4			
5	Chair Force One Q2		
6	Finalize Requirements & Refinement	05-Jan-2026	11-Jan-2026
7	Review Quarter 1 Design	08-Jan-2026	08-Jan-2026
8	Identify Critical Test to Evaluate Functional Performance	08-Jan-2026	08-Jan-2026
9	FEA Validation & Analytical Checks of Quarter 1 Design	09-Jan-2026	10-Jan-2026
10	Project Risk Review to Identify Technical, Safety, and Budget Challenges	11-Jan-2026	11-Jan-2026
11	Manufacturing Preparation	13-Jan-2026	02-Feb-2026
12	Finalize Sub Systems - Attachments, L.A, Controller, and Wheels	13-Jan-2026	13-Jan-2026
13	Create Manufacturing Plans - Built to Order vs PRL	15-Jan-2026	29-Jan-2026
14	Solidify CAD File and Test Sub System Parts	19-Jan-2026	02-Feb-2026
15	Functional Prototype	20-Jan-2026	02-Feb-2026
16	3D Printing Pieces that May Need More Testing Before Manufacturing	20-Jan-2026	25-Jan-2026
17	Send Parts out for Manufacturing and/or Build in the PRL	20-Jan-2026	02-Feb-2026
18	Actuator Integration - Test Synchronization and Electronic Components	27-Jan-2026	30-Jan-2026
19	Build Testing Fixtures to Test the Full Scale Prototype	30-Jan-2026	02-Feb-2026
20	Testing Functional Prototype	02-Feb-2026	22-Feb-2026
21	Engineering Requirements Testing - Based on Final ER List	03-Feb-2026	10-Feb-2026
22	Redesign Weak Points After Testing	08-Feb-2026	22-Feb-2026
23	Update Manufacturing Plans After Resolving Weak Points	10-Feb-2026	22-Feb-2026
24	Final Prototype	16-Feb-2026	28-Feb-2026
25	Manufacture Final Parts - Build to Order and/or PRL	16-Feb-2026	28-Feb-2026
26	Assemble Prototype to Ensure It's Ready for the Final Presentation	23-Feb-2026	28-Feb-2026
27	Final Evaluations and Presentation	28-Feb-2026	12-Mar-2026
28	Final Documentation - All Components	02-Mar-2026	10-Mar-2026
29	Final Report	07-Mar-2026	10-Mar-2026
	Final Presentation	10-Mar-2026	12-Mar-2026
~~			



Q2 Budget Overview

- Spent: \$362
- **Anticipated Total Spend: \$2240**
- Major Anticipated Expenses
 - Linear actuators: ~\$720
 - 2x foldable wheelchair + scissor brakes: ~\$500
 - Stock/hardware for parts + test fixtures: ~\$515



Conclusion

- **Objective:** Design a wheelchair that allows manual wheelchair users to independently board and exit an aircraft, move through its aisle, and seat themselves on standard aircraft seating.
- **Recommended Solution:** A manual, foldable, dual-mode wheelchair that stows in an overhead bin.



Alfonse the Duck
Fully Approves this Concept

Backup Slides

Selection Rationale

- **Our Users**
 - What's important to them
- **Engineering Requirements**
 - Feasibility
- **Potential Risks**
 - Stowing
 - Modifications


Modular Attachment vs Full Solution

Applicable ERs	Weight	Modular Attachment	Full Solution
ER 1-5: Airplane Stowability	2	1	3
ER 1-1: End-User Independence	3	1	2
ER 3-3a, b: Durability	1	1	3
ER 4-2a, b, c: Safety to End-User	3	2	2
	Total	12	21



Decision matrices 3 - Rigid vs. Foldable Wheelchair

Applicable ERs	Priority (1-3 scale)	Rigid	Foldable
Can be independently operated by the user ER 1-1a	3	3	3
Without our attachment - wheelchair can brake ER 1-1b	1	3	2
Fits through the cabin aisle ER 1-2	3	3	3
Can get over the gap between the plane and jet bridge ER1-3	2	2	2
Can be stowed away ER 1-5	3	1	3
Must support 300 lbs ER 3-1a	3	3	3
Must weight under 50 lbs ER 4-1d	3	3	3
	Total	46	51

The Linear Actuator Circuit






Attachment of Linear Actuator to Wheelchair




Wheelchair Folded

- Required volume to fit in: $54 \times 22 \times 14$ inches
- Wheelchair w/o parts: $\sim 27 \times 21 \times 12$ inches
- Main wheels are 24 inch diameter, can fit along diagonal of volume (~ 26 inch)
 - Can also reduce diameter using other COTS wheels



Proposed Stowing Configuration



Design/Stress Analysis of Backrest

Main risk: Tube stress concentrations at edges



Tube Connector

$$M = Fd = (52.5 \text{ lb}f)(7 \text{ in}) = 367.5 \text{ lb-in}$$

$$I = \frac{\pi}{64} d_o^4 = 0.01655 \text{ in}^4$$

$$\sigma_{max} = \frac{Md_o}{2I} \approx 8460 \text{ psi}$$

$$\sigma_{yield, 6061} \approx 40,000 \text{ psi} > \sigma_{max}$$



ASME Canons Incorporated into Design

- (1) Engineers shall **hold paramount the safety, health, and welfare of the public** in the performance of their professional duties.

- (4) Engineers shall **act in professional matters for each employer or client as faithful agents or trustees** and shall avoid conflicts of interest or the appearance of conflicts of interest.

Ethics of the Design

- Considerations
 - Safety for primary and secondary users
 - Balancing interests of Boeing, airlines, and wheelchair users
 - Comfort & independence
 - Efficiency (time & cost)
- Steps Taken
 - Design to sufficient factor of safeties and multiple failure modes
 - Wheelchair max weight, 50 lbs
 - Improved independence of end-user
 - Linear actuators
 - Stowable in an overhead bin
 - Reduces boarding/deboarding times